
HL7 Conformance Profile

Oru Mediweb

Copyright notice:

Copyright 2024 AGFA HealthCare
All rights reserved

Agfa, the Agfa rhombus, Point of Knowledge, and See More. Do More, ... (other trademarks) are trademarks of Agfa Gevaert N.V., Belgium or its affiliates. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringement.

The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications which must be met by Agfa. All information contained herein is intended for guidance purposes only, and characteristics of the products described in this publication can be changed at any time without notice.

Products may not be available for your local area. Please contact your local sales representative for availability information.

Agfa diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

Publication date:

29 03 2024

Publication number:

Not yet assigned

Corporate address:

AGFA HealthCare
SEPTESTRAAT 27
B-2640 MORTSEL
BELGIUM
+32(3)4448400

Conformance parameters

Message Profile

- HL7 Version: 2.4
- Profile Type: HL7
- Topics: confsig-Agfa-2.4-profile-accNE_accAL-Deferred

Encoding Method

ER7

Interaction 1

Dynamic Definition

- Accept Acknowledgement: NE
- Application Acknowledgement: NE
- Acknowledgement Mode: Immediate

Static Definition

- Message Type: ORU
- Trigger Event: R01
- Message Structure: ORU_R01
- Topics: confsig-Agfa-2.4-static-ORU-R01-null-ORU_R01-2005.1.4--Sender

Message structure

MSH { [PID [[PV1]]] { [ORC] OBR { [NTE] } { [OBX] { [NTE] } } } }

MSH - Message Header

- Usage: Required
- Cardinality:1..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
1	Field Separator	ST		1	R	1..1	
2	Encoding Characters	ST		4	R	1..1	
3	Sending Application	HD	HL70361	227	O	0..1	
3.1	namespace ID	IS	HL70363	180	O	..	
3.3	universal ID type	ID	HL70301	3	C	..	
4	Sending Facility	HD	HL70362	227	O	0..1	
4.1	namespace ID	IS	HL70363	180	O	..	
5	Receiving Application	HD	HL70361	227	O	0..1	
5.1	namespace ID	IS	HL70363	180	O	..	
6	Receiving Facility	HD	HL70362	227	O	0..1	
6.1	namespace ID	IS	HL70363	180	O	..	
7	Date/Time Of Message	TS		26	R	1..1	
7.1	Date/Time	NM		14	R	..	e.g. 20040328134602
8	Security	ST		40	O	0..1	
9	Message Type	CM_MSG	HL70076	15	R	1..1	
9.1	message type	ID	HL70076	3	O	..	
9.2	trigger event	ID	HL70003	3	O	..	
9.3	message structure	ID	HL70354	7	O	..	
10	Message Control ID	ST		20	R	1..1	
11	Processing ID	PT		3	R	1..1	
11.1	processing ID	ID	HL70103	3	O	..	
12	Version ID	VID	HL70104	973	R	1..1	
12.1	version ID	ID	HL70104	3	O	..	

3.1. namespace ID

parameter: MSH_SendingApplic, value: 'QDOC'

4.1. namespace ID

parameter: MSH_SendingFacility, value: 'AGFA'

5.1. namespace ID

parameter: MSH_ReceivingApplic, value: 'RAPP'

6.1. namespace ID

parameter: MSH_ReceivingFacility, value: 'RFAC'

7.1. Date/Time

YYMMDDHHMMSS

Current date

9.1. message type

fixed value: 'ORM'

9.2. trigger event

fixed value: '001'

9.3. message structure

fixed value: 'ORM_001'

Segment group: PATIENT RESULT

- Usage: Required

- Cardinality: 1..*

Segment group: PATIENT

- Usage: Optional

- Cardinality: 0..1

PID - Patient identification

- Usage: Required

- Cardinality: 1..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
1	Set ID - PID	SI		20	O	0..1	
2	Patient ID	CX		1913	O	0..1	
2.1	ID	ST		20	O	..	
3	Patient Identifier List	CX		1913	R	1..*	
3.1	ID	ST		200	R	..	
3.4	assigning authority	HD		407	O	..	
3.4.1	namespace ID	IS	HL70363	200	O	..	
3.5	identifier type code (ID)	ID	HL70203	200	O	..	
3.6	assigning facility	HD		208	O	..	
3.6.1	namespace ID	IS	HL70363	200	O	..	
3.6.2	universal ID	ST		3	O	..	
3.6.3	universal ID type	ID	HL70301	3	O	..	
4	Alternate Patient ID - PID	CX		1913	O	0..1	
4.1	ID	ST		20	O	..	
5	Patient Name	XPN		1103	O	0..*	
5.1	family name	FN		46	O	..	
5.1.1	surname	ST		30	O	..	
5.1.2	own surname prefix	ST		3	O	..	

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
5.1.3	own surname	ST		3	O	..	
5.1.4	surname prefix from partner/spouse	ST		3	O	..	
5.1.5	surname from partner/spouse	ST		3	O	..	
5.2	given name	ST		30	O	..	
5.3	second and further given names or initials thereof	ST		30	O	..	
5.4	suffix (e.g., JR or III)	ST		3	O	..	
5.5	prefix (e.g., DR)	ST		3	O	..	
6	Mother's Maiden Name	XPN		1103	O	0..*	
6.1	family name	FN		46	O	..	
6.1.1	surname	ST		30	O	..	
6.1.2	own surname prefix	ST		3	O	..	
6.1.3	own surname	ST		3	O	..	
6.1.4	surname prefix from partner/spouse	ST		3	O	..	
6.1.5	surname from partner/spouse	ST		3	O	..	
6.2	given name	ST		30	O	..	
6.7	name type code	ID	HL70200	1	O	..	
7	Date/Time Of Birth	TS		26	O	0..1	
7.1	Date/Time	NM		14	R	..	e.g. 20040328134602
8	Administrative Sex	IS	HL70001	3	O	0..1	
11	Patient Address	XAD		631	O	0..*	
11.1	street address (SAD)	SAD		328	O	..	
11.1.1	street or mailing address	ST		264	O	..	
11.3	city	ST		30	O	..	
11.4	state or province	ST		40	O	..	
11.5	zip or postal code	ST		20	O	..	
11.6	country	ID	HL70399	40	O	..	
13	Phone Number - Home	XTN		850	O	0..*	
13.7	Phone number	NM		31	O	..	
18	Patient Account Number	CX		1913	O	0..1	
18.1	ID	ST		50	O	..	
18.4	assigning authority	HD		387	O	..	
18.4.1	namespace ID	IS	HL70363	180	O	..	
19	SSN Number - Patient	ST		50	O	0..1	
23	Birth Place	ST		30	O	0..1	
29	Patient Death Date and Time	TS		26	O	0..1	
29.1	Date/Time	NM		14	R	..	e.g. 20040328134602

2.1. ID

patients.p_extracode, patients.p_code

3.1. ID

patient_identifiers.pid_id, patients.p_natnumber

3.4.1. namespace ID

patient_id_assigning_auth.pidasauth_namespaceid

3.5. identifier type code (ID)

patient_id_assigning_auth.pidasauth_namespaceid

3.6.1. namespace ID

patient_id_assigning_facility.pidasfac_namespaceid

4.1. ID

patients.p_extracode

5.1.1. surname

patients.p_name

5.2. given name

patients.p_fname

5.3. second and further given names or initials thereof

patients.p_mname

6.1.1. surname

patients.p_name

6.2. given name

patients.p_fname

6.7. name type code

fixed value: 'B'

7.1. Date/Time

YYYYMMDDHHMMSS

patients.p_born

8. Administrative Sex

patients.p_sex

11.1.1. street or mailing address

patients.p_street & patients.p_streetnr

11.3. city

patients.p_town

11.4. state or province

patient_items['PROVINCE'].pi_value

11.5. zip or postal code

patients.p_zip

11.6. country

patients.p_country

13.7. Phone number

patients.p_phone ~ patients.p_phone

19. SSN Number - Patient

patients.p_siscode

23. Birth Place

patients.p_birthplace

29.1. Date/Time

YYYYMMDDHHMMSS

patients.p_datedied

Segment group: VISIT

- Usage: Optional
- Cardinality: 0..1

PV1 - Patient visit

- Usage: Optional
- Cardinality: 0..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
19	Visit Number	CX		1913	O	0..1	
19.1	ID	ST		20	O	..	

19.1. ID

used in select service_episodes.se_admission_number

End of segment group VISIT

End of segment group PATIENT

Segment group: ORDER OBSERVATION

- Usage: Required
- Cardinality: 1..*

ORC - Common Order

- Usage: Optional
- Cardinality: 0..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
13	Enterer's Location	PL		1230	O	0..1	
13.1	point of care	IS	HL70302	100	O	..	

13.1. point of care

services.s_hospitalkey

OBR - Observation Request

- Usage: Required
- Cardinality: 1..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
1	Set ID - OBR	SI		3	O	0..1	
3	Filler Order Number	EI		427	O	0..1	
3.1	entity identifier	ST		64	O	..	
3.2	namespace ID	IS	HL70363	64	O	..	
4	Universal Service Identifier	CE		483	O	0..1	
4.1	identifier	ST		3	O	..	
4.2	text	ST		254	O	..	
4.3	name of coding system	IS	HL70396	64	O	..	
4.4	alternate identifier	ST		16	O	..	
7	Observation Date/Time #	TS		26	O	0..1	
7.1	Date/Time	NM		14	O	..	e.g. 20040328134602

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
13	Relevant Clinical Info.	ST		254	O	0..1	
14	Specimen Received Date/Time *	TS		26	O	0..1	
14.1	Date/Time	NM		14	O	..	e.g. 20040328134602
15	Specimen Source	CM_SPS	HL70070	4436	O	0..1	
15.1	specimen source name or code	CE		717	O	..	
15.1.1	identifier	ST		254	O	..	
16	Ordering Provider	XCN		3002	O	0..1	
16.1	ID number (ST)	ST		20	O	..	
16.2	family name	FN		244	O	..	
16.2.1	surname	ST		100	O	..	
16.3	given name	ST		30	O	..	
22	Results Rpt/Status Chng - Date/Time +	TS		26	O	0..1	
22.1	Date/Time	NM		14	O	..	e.g. 20040328134602
24	Diagnostic Serv Sect ID	ID	HL70074	254	O	0..1	
25	Result Status +	ID	HL70123	10	O	0..1	
26	Parent Result +	CM_PRL		755	O	0..1	
26.1	OBX-3 observation identifier of parent result	CE		11	O	..	
26.1.1	identifier	ST		3	O	..	
26.1.2	text	ST		3	O	..	
27	Quantity/Timing	TQ		1545	O	0..1	
27.6	priority	ST		3	O	..	
32	Principal Result Interpreter +	CM_NDL		835	O	0..1	
32.1	name	CN		147	O	..	
32.1.1	ID number (ST)	ST		10	O	..	
33	Assistant Result Interpreter +	CM_NDL		835	O	0..1	
33.1	name	CN		565	O	..	
33.1.2	family name	FN		254	O	..	
33.1.3	given name	ST		254	O	..	

3.1. entity identifier

lab_service_request.lsr_service_request_id

3.2. namespace ID

lab_service_request.lsr_filler_id

4.2. text

lab_topic.lto_name

4.3. name of coding system

lab_topic.codesys

4.4. alternate identifier

lab_topic.lto_position

7.1. Date/Time

YYYYMMDDHHMMSS

lab_sample.ls_collectiondate

13. Relevant Clinical Info.

lab_service_request_ci.lsrc_value,
lab_blob.lb_value

14.1. Date/Time

YYYYMMDDHHMMSS

lab_sample.ls_sampleid, lab_sample.ls_receivedate

15.1.1. identifier

lab_sample.ls_receivedate, lab_sample_ci.lsc_value

16.1. ID number (ST)

doctors.d_hospitalkey

16.2.1. surname

doctors.d_lname

16.3. given name

doctors.d_fname

22.1. Date/Time

YYYYMMDDHHMMSS

lab_service_request.lsr_versiondate,

lab_sample.ls_versiondate

24. Diagnostic Serv Sect ID

lab_service_request_ci.lsrc_value

25. Result Status +

lab_sample.ls_status,

lab_service_request.lsr_status

32.1.1. ID number (ST)

lab_sample.ls_resp_doctor

33.1.3. given name

lab_sample_ci.lsc_value

NTE - Notes and Comments

- Usage: Optional

- Cardinality:0..*

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
3	Comment	FT		254	O	0..1	
4	Comment Type	CE	HL70364	483	O	0..1	
4.1	identifier	ST		20	O	..	

3. Comment

lab_sample_ci.lsc_value, arch_lab_blob.alb_value

4.1. identifier

lab_class_item_def.lci_code

Segment group: OBSERVATION

- Usage: Required

- Cardinality:1..*

OBX - Observation/Result

- Usage: Optional

- Cardinality: 0..1

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
2	Value Type	ID	HL70125	2	O	0..1	
3	Observation Identifier	CE		483	O	0..1	
3.1	identifier	ST		64	O	..	
3.2	text	ST		254	O	..	
3.3	name of coding system	IS	HL70396	64	O	..	
3.4	alternate identifier	ST		16	O	..	
4	Observation Sub-Id	ST		10	O	0..1	
5	Observation Value	*		3	O	0..1	
6	Units	CE		483	O	0..1	
6.1	identifier	ST		254	O	..	
7	References Range	ST		254	O	0..1	
8	Abnormal Flags	IS	HL70078	10	O	0..1	
11	Observation Result Status	ID	HL70085	10	O	0..1	
14	Date/Time of the Observation	TS		26	O	0..1	
14.1	Date/Time	NM		14	O	..	e.g. 20040328134602

2. Value Type

fixed values: 'DT', 'FT', 'CE', 'ST'

3.1. identifier

lab_topic.lto_code

3.2. text

lab_topic.lto_name

3.3. name of coding system

lab_topic.lto_codesys

3.4. alternate identifier

lab_topic.lto_position

4. Observation Sub-Id

lab_result.lr_test_sub_id

6.1. identifier

lab_measurement_def.lm_value

7. References Range

lab_measurement_def.lm_value

8. Abnormal Flags

fixed values: 'R', 'S', 'T', 'VS', 'MS' lab_result.lr_abnormal

11. Observation Result Status

fixed values: 'T', 'F', 'D', 'X' lab_result.lr_status

14.1. Date/Time

YYYYMMDDHHMMSS

lab_result.lr_date

NTE - Notes and Comments

- Usage: Optional
- Cardinality:0..*

Seq.	Name	Type	Table	Len.	Opt.	Card.	Contents
3	Comment	FT		3	O	0..1	

3. Comment

lab_blob.lb_value

End of segment group OBSERVATION

End of segment group ORDER_OBSERVATION

End of segment group PATIENT_RESULT